


MEDIWISS Analytic GmbH

Uerdinger Strasse 3
D-47441 Moers
02841/8890470
pwahl@mediwiss-analytic.de



AllergyScreen™/AlleisaScreen™
CIRCULAR LETTRE 8

RDFs to use for your connect cameras-
(depending on the software-version and the date of camera-production)

	Connects - Produced before May 2009	Connects - Produced in May 2009 or later
Software 490	v13	---
Software 510	v13	---
Software 520	v13	---
Software 530	v13	---
Software 540	v23	v40

Important information:
The date of camera- production is mentioned in the serial-number.

SN **0509**1000xxxx means: produced in May 2009